



The Clatterbridge Cancer Centre

NHS Foundation Trust



University Hospitals
Coventry and Warwickshire

NHS Trust



Background Neutron Scatter Modelling in Proton Beam Therapy

Sarah Baker


16th June 2016

Talk Overview

- Project Rationale and Outline
- Equipment
- Results
 - Part 1: Neutron response of a plastic scintillator
 - Part 2: Comparison to neutrons from californium 252
 - Part 3: Neutron production at Clatterbridge Cancer Centre
- Conclusions and Future Plans

Rationale: Proton Therapy

- Advanced radiotherapy using protons, rather than X-rays
- Maximum dose in the final few millimetres of the proton tracks (Bragg peak)
 - Bragg peak position dependent on the energy of the incident proton
 - Sharp distal dose fall-off
- Advantages:
 - Treatments close to critical organs such as head and necks
 - Paediatrics - reduction of induced secondary cancers in later life



Rationale: Neutron Production

- Proton beams interact with the beam-line and within the patient to produce secondary radiation
 - Neutrons can cause significant biological damage
- Neutron damage to tissues is strongly dependent on neutron energy and flux
- Spectroscopy of secondary neutrons would permit reduction of neutron scatter dose
 - Optimisation of beam-line arrangements
 - Improvements to dose delivery methods
 - Shielding of sensitive equipment
 - Reduction of risk of neutron induced secondary cancers

Project Outline: Remit

- UCL have developed a scintillator to measure the energy resolution and spread of a proton beam:

An optical module designed for the SuperNEMO experiment has measured the 60 MeV proton beam at Clatterbridge with an energy resolution of:


1.58 ± 0.27 % FWHM
(for 40 MeV “visible” energy)

simple extrapolation:

0.7 % FWHM @ proton imaging energy (300 MeV)

Limit on energy spread of Clatterbridge beam:

$\text{FWHM (60 MeV)} < 1.47$ % @ 90% CL



Source: Seminar by Anastasia Basharina-Freshville, UCL

Project Aim: Model the scintillator response to neutrons


Project Outline: Particle Interactions

- Neutrons uncharged so interact with atomic nuclei
 - Scattering (elastic and inelastic)
 - Absorption (fission and capture)
- Neutron cross sections are strongly dependent on energy
 - Typically higher at lower energies
 - Detection often requires 'moderators' to reduce neutron energy via scattering – hydrogen is particularly effective
- Full simulation requires other relevant particle interactions (proton, γ , e^- , e^+ ...)
- Detector scintillation properties and resolution effects must also be considered

Equipment: Monte Carlo Simulation


- Simulation of physical processes using random sampling of probability distributions
- Monte Carlo in medical physics:
 - Event generation
 - Beam-line and detector simulation
- Geant4 used to generate input particle spectra and to model the scintillator
 - Applications in particle, nuclear, accelerator and dark matter physics
 - Provides access to particle data libraries
 - Must consider which process models and energies are important

Equipment: Scintillator and PMT




Wrapping:
Sides: 75 μm of PTFE (Teflon) ribbon
Sides and entrance face: 12 μm of Mylar

- Plastic scintillator Polyvinyl toluene (PTV) hexagonal block
- Wrapped in reflective material
- Coupled to a photo-multiplier tube (PMT)
- Irradiated with beam of particles




Source: Seminar by Anastasia Basharina-Freshville, UCL

Results Part 1: Monoenergetic Beam




- Peak at 1 MeV - neutron energy fully deposited in scintillator (hydrogen scattering)
- Scatter from C12 atom
- 2.2 MeV gamma produced in neutron capture for hydrogen

Results Part 1: Polyenergetic Beam




- Peak at 1 MeV - neutron energy fully deposited in scintillator (hydrogen scattering)
- Scatter from C12 atom
- 2.2 MeV gamma produced in neutron capture for hydrogen

Results Part 1: Polyenergetic Beam



- Peak at 1 MeV - neutron energy fully deposited in scintillator (hydrogen scattering)
- Scatter from C12 atom
- 2.2 MeV gamma produced in neutron capture for hydrogen

Results Part 1: Quenched Energy




- Light yield is non-linear (quenched) - Birks law:

$$\frac{dY}{dx} = \frac{S}{1 + kB(dE/dx)} \times \frac{dE}{dx}$$

dY/dx light yield per unit path length
 dE/dx energy lost by particle per unit path length
 kB relates density of ionisation to energy loss
 S absolute scintillation efficiency


- Neutrons transfer energy to C12 and hydrogen

Results Part 2: CF-252 Gamma/Neutron Multiplicity (FREYA)




- Test Monte Carlo using a Californium source
- CF-252 emits gammas and neutrons via spontaneous fission
- Spontaneous fission not well handled by Geant4
- Alternative generator data (FREYA) used to create a bespoke CF-252 generator

Results Part 2: CF-252 Gamma Spectrum (FREYA)





- Multiplicities and spectra combined to create a probability based particle generator to recreate the spontaneous fission
- Does not account for correlations

Results Part 2: CF-252 Neutron Spectrum (FREYA)




- Multiplicities and spectra combined to create a probability based particle generator to recreate the spontaneous fission
- Does not account for correlations


Results Part 2: CF-252 Generator Spectrum



Results Part 2: CF-252 Deposited Energy




Results Part 2: CF-252 Quenched Energy




- Peak at neutron ~ 2 MeV in CF-252 visible energy spectrum
- Effect of the quench modelling?

Results Part 2: CF-252 Deposited Energy




- Peak at neutron ~ 2 MeV in CF-252 visible energy spectrum
- Effect of the quench modelling?
- Quenching of the 2.2 MeV gamma from neutron capture?

Results Part 2: CF-252 Deposited Energy




- Monte Carlo – therefore can remove neutron capture processes from the simulation entirely (not physical)
- No peak seen

Results Part 2: CF-252 Data




- Measured by UCL using scintillator
- Peak also seen at 2 MeV - as seen in MC?

Results Part 2: CF-252 Data/MC Comparison




- Total energy deposited
- Shape match only (data black, MC grey)
- Deficit in MC events below 2 MeV
- Excess in MC events at lowest energies

Results Part 2: Threshold Effects at Low Energy




- Detector has a minimum energy threshold for detection
- Causes a turn on curve that depends on the value of the threshold
- Statistical errors too small to view (8 MeV bump not due to statistics)

Results Part 2: Model Validation



- Lack of MC events below 2 MeV could be due to the generator modelling
- Comparison to Lawrence Livermore National Laboratory (LLNL) fission package
- Agreement in neutrons, disagreement in gammas

Results Part 2: Improve the Neutron Source




- Lead preferentially attenuates gammas
- Measurements with lead shielding between CF-252 and scintillator will mean a purer neutron source
- Reduce impact of uncertainty of gamma modelling

Beam Line at Clatterbridge




- Ultimate goal – develop a neutron spectrometer
- Provide in situ measurements of neutron dose
- Optimise beam-line to reduce dose to patients and sensitive equipment

Results Part 3: Modelling The Beam Line in Geant4




- Core beam-line components modelled
- Concrete walls, ceiling and floor also added to recreate scatter conditions (not shown)
- 62.5 MeV protons (reduced to ~60 MeV by air scatter)
- 1,000,000 events

Results Part 3: Production of Neutrons in the Room





- Model (currently) shows about 7 neutrons produced per 100 protons in the room.
- As model complexity increases, more things to interact with will result in more neutrons
- These neutrons scatter widely

Results Part 3: Energy deposition in the Clatterbridge Treatment room




Results Part 3: Energy deposition in the Clatterbridge Treatment room



Source: Howell, Med Phys. 2014 Sep; 41(9): 092104


250 MeV Beam at 100cm from Isocentre.

Results Part 3: Internal Neutron Production



- Neutrons are produced within the patient and detector and contribute to neutron dose ~ 1 neutron per 100 protons
- A correction may be required if patient dose calculations are performed using scintillator data

Results Part 3: Internal Neutron Production



- Neutrons are produced within the patient and detector and contribute to neutron dose ~ 1 neutron per 100 protons
- A correction may be required if patient dose calculations are performed using scintillator data

Summary and Future Work

- The response of the detector to neutrons has been modelled using Geant4 and compared to a CF-252 source
- The beam-line in Clatterbridge has been modelled and the neutron production within the room assessed
- Future plans:
 - CF-252 comparison with lead shielding
 - Dose deposition as a function of distance from beam-line
 - Perform test-beam measurements at Clatterbridge to validate beam-line model and detector simulation



Thanks for Listening!



Clatterbridge Beam-Line