

Proton Beam Therapy


HEP Group Day 2025

Sonia Escribano-Rodriguez, Joe Bateman, Raffaella Radogna, Febian Febian, Simon Jolly



Proton beam therapy

Proton Beam Therapy (PBT) offers a highly localised depth-dose distribution, due to the well-defined range of protons in matter.






Quality Assurance (QA)

To ensure accurate dose delivery in PBT, Quality Assurance (QA) measurements are required.

- Current methods often compromise between speed and accuracy.
- Measurements of <u>beam range</u>, absolute dose, spot position and size.

Develop a detector to make fast, accurate measurements of proton range.






QuARC

QuARC = **Quality Assurance Range Calorimeter**

- Optically isolated polystyrene scintillator sheets of 100 mm x 100 mm x 3 mm.
- Each sheet is coupled to a photodiode. PD read-out by an FPGA.
- Analytical model deployed to fit depth-light curve to reconstruct Bragg curve.







What's new in PBT?


- New polished scintillators
- Photodiodes on both sides of the scintillator
- Improved mechanical set-up









 Testing response of the new polished scintillator sheets with shoot-through measurements.








 Testing response of the new polished scintillator sheets with shoot-through measurements.









Focusing on module 2, made of polished sheets we also see a different light response for ST front and ST back.


When the board is placed on the left it always sees more light -> beam misalignment?









Looking at the actual BPs measured:








Confirm if there is a beam misalignment at PPTC or if it's a detector issue






Shift the beam with respect to the centre of the detector to study misalignment issues from Prague







Shift the beam with respect to the centre of the detector to study misalignment issues


$$shift[i] = \frac{PD_A[i] - PD_B[i]}{PD_A[i] + PD_B[i]}$$
 PPTC beam is shifted ~1cm to the left













Next Steps

QuARC → QuADPro

Quality Assurance Range Calorimeter → Quality Assurance Detector for Protontherapy

- Integrate a fibre array to measure beam position and spot size
- Beam test in TIFPA @FLASH rates next week!
- Add photodiode boards to 4 sides of the scintillator
- Upgrade to the electronics







+ Harry Barnett and Matt Warren

